



PMI SCIENCE
PHILIP MORRIS INTERNATIONAL

Lung/liver-on-a-chip system

Dr. David Bovard, Scientist - Microphysiological Systems

**PMI R&D, Philip Morris Products S.A.,
Quai Jeanrenaud 5, CH-2000 Neuchâtel, Switzerland**

Philip Morris International R&D facilities

Neuchatel, Switzerland




Singapore



- 200+ world-class scientists, engineers, and medical directors (40+ nationalities)
- 30+ disciplines, including material science, consumer electronics, *in vitro* and *in vivo* research, computer science, and clinical science

Aerosol absorption



Lung

- CYP 1A1/1B1
- CYP 2B6
- CYP 2E1
- CYP 2J2
- CYP 3A5

Liver

- CYP 1A2
- CYP 2E1
- CYP 2C8
- CYP 2C9
- CYP 3A4


Key features of the platform




- Air-liquid interface (ALI) bronchial model
- Liver spheroid model

Lung and liver models


Bronchial ALI model



NHBE cells



Liver spheroid model




HepaRG™ cells





Model stability

- Key tissue characteristics were assessed over a six-week period


Bronchial ALI model



Liver spheroid model



Key features of the platform




Lung/liver-on-a-chip platform

Pump controller




Pump main unit




Chip

Reservoir plate


Polyetheretherketone:
Biocompatible and
nonabsorbent material



After eight hours of incubation, nicotine concentration remained stable in the PEEK chip but not in the PDMS chip.




Bronchial ALI model




Liver spheroid model

Key features of the platform




Defined coculture medium

Bronchial ALI model




Liver spheroid model




10 days



William's E Medium




PneumaCult™-ALI Medium




	24 h	48 h	72 h	96 h	168 h	240 h
CYP3A4	18.1 ***	47.0 ***	19.7 **	53.1 ***	31.5 ***	18.2 ***
CYP3A43	6.9 *	33.8 *	9.3 *	8.5 ***	2.8	5.5 *
CYP2B6	3.9 *	4.8 *	2.0	4.6 **	4.3 *	16.7 **
CYP7A1	3.2 **	3.3 **	1.7 *	2.5 **	1.1	1.1
CYP2C8	3.1 **	1.6	2.1 ***	2.3 *	2.7 ***	1.1
CYP2A13	3.0 *	1.9 *	1.7	8.5 *	6.0 **	3.8 *
ADH7	2.4	-1.4	-1.9	-1.4 *	-2.2	-6.8
CYP7B1	1.9 **	6.1 **	1.2	3.3 ***	2.6 **	-1.3
CYP2C9	1.8 *	2.9 **	1.3	4.3 **	2.6 ***	-4.2
ADH1A	1.8	1.4	1.1	2.7 *	1.3	-1.5
CYP2F1	1.8 *	1.1	-1.4	1.9	-1.0	1.4
CYP4F12	-1.8 *	-2.1 **	-2.3 **	1.3	-1.7 **	1.1
AADAC	-1.8 **	1.3	-2.2 **	1.2	-1.6	-2.7 **
ALDH5A1	-1.9	-2.7 **	-6.5	1.5	-1.7	2.5 *
CYP27A1	-1.9 *	-1.1 *	-2.0 *	1.8 **	1.1	-1.1
CYP4A22	-2.0	-3.3 *	-3.6	-1.1	-1.4	-1.7 *
CYP1B1	-2.1 *	-2.0 *	-2.7 **	1.3	-1.3	-1.0
CYP27B1	-2.3	1.8	-1.8	5.2 *	1.1	-1.7
CYP26B1	-2.8 **	-2.5 *	-5.8 **	1.0	-1.0	1.2
CYP1A1	-3.2 ***	-3.6 ***	-5.0 ***	-2.7 *	-1.4	1.8
CYP19A1	-3.7 *	-1.9	-3.0 *	-1.9	-3.3 **	-2.1 *
ALDH3A1	-4.1	-4.6 **	-7.0 *	1.1	1.1	1.1
CYP4A11	-4.9 ***	-4.8 ***	-9.0 **	-1.8 *	-4.9 ***	-3.2 **
ALDH1A3	-10.1 *	1.1	-5.9	1.3	-1.6	1.2






Key features of the platform




Lung-liver coculture






- Tissues were maintained for 28 days in either
 - Static (monoculture, in standard 24-well plate)
 - Dynamic (monoculture, in the chip with medium flow)
 - Coculture (coculture with liver spheroids)




Lung-liver coculture



- Tissues were maintained for 28 days in either
 - Static (monoculture, in standard 96-well plate)
 - Dynamic (monoculture, in the chip with medium flow)
 - Coculture (coculture with bronchial ALI tissues)




Key features of the platform




AFB1

- Produced by *Aspergillus flavus*
- Common food contaminant
- Causes hepatocellular carcinoma



Liver cancer data from the GLOBOCAN 2002 database
(http://www-dep.iarc.fr/GLOBOCAN_frame.htm)

Aflatoxin data from Williams et al., *Human Aflatoxicosis in Developing Countries*,
Am J Clin Nutr 80:1106–22, 2004.






AFB1 in the chip

- AFB1: 100 μM
- Compound diluted in medium
- 48-hour exposure



AFB1 in the chip

- Results:
 - AFB1 caused strong cytotoxicity to the bronchial ALI tissues
 - AFB1 cytotoxicity to the ALI bronchial cultures was decreased by the presence of liver spheroids in the chip



Conclusion & perspectives

- A two-organ chip has been developed using a biocompatible and nonabsorbent material
- Lung and liver tissues can be maintained for 28 days in coculture in the chip
- The lung-liver crosstalk was demonstrated using AFB1
- **The lung-liver-on-a-chip platform allows for chronic and acute toxicity testing**
- Using the current platform, several aspects will be optimized:
 - Higher throughput
 - Increased tissue complexity (cell type composition)
 - Connection of additional “organs”

Thank you for your attention !!

Stefan
Frentzel

Moran
Morelli

Arno
Knorr

Keyur
Trivedi

Julia
Hoeng

Diego
Marescotti

Anita
Iskandar

Antonin
Sandoz

Emmanuel
Guedj

Karsta
Luettich